

Try Our New Sound System! Hearing Aid Wearers Will Not Believe Their Ears!



HEARING LOOP INSTALLED
Switch hearing aid to T-coil

Player's Theatre has added a Hearing Loop, which is the only system that sends the voices of the performers from our microphones directly into a hearing aid or cochlear implant. The system provides the clearest, cleanest sound imaginable and does not require a headset, as it works with the **T-Coil** or Telecoil that most hearing aids come equipped with. **The system is available to anyone sitting in Rows A to K.**

A T-coil is a tiny wireless receiver inside hearing aids and it is accessible via the programming button or switch on your aid. If you have never used your T-Coil, it is recommended that you consult your audiologist to confirm it has been activated and to learn how to access it. If you cannot access your T-coil or do not have hearing aids, please ask the usher for a Loop Receiver which you can use to receive the system at today's event. Again, the Loop Receiver is only for use in rows A to K. You will be amazed at how much more you can hear and understand, given the clarity of sound!

If you have more questions about hearing loops or T-coils, please contact the Hearing Loss Association of Sarasota at 941-706-4312.